

# Sondrelix — Stability Protocol

Protocol STAB-P-2024-10

## Sample Inventory and Pull Windows

Sample inventory per condition and timepoint is defined in the pull schedule appendix. Scheduled pulls are performed within  $\pm 7$  days of the nominal date; pull records are maintained in the site stability management system.

## Container Closure and Orientation

Samples are stored in the commercial container closure system. Bottles are stored upright; a subset of inverted samples is included at the initial and final timepoints to assess closure-liner interaction. Induction seal integrity is verified at each pull.

## Storage Conditions

The primary (intermediate/zone IVa condition) storage condition for this protocol is  $30^{\circ}\text{C} \pm 2^{\circ}\text{C}$  /  $65\% \text{ RH} \pm 5\% \text{ RH}$ . Accelerated studies at  $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$  /  $75\% \text{ RH} \pm 5\% \text{ RH}$  were completed through six months during development and are not part of this protocol.

## Pull Schedule

Long-term samples are tested at 0, 3, 6, 9, 12, 18, 24 and 36 months.

## Post-Approval Commitments

As a condition of approval, three consecutive production batches of Sondrelix will be entered into the stability program, and the retest period will be re-evaluated after 36-month data are available.

## Bracketing Justification

Bracketing is justified under ICH Q1D on the basis that the 25 mg and 100 mg presentations bracket the 75 mg presentation with respect to fill volume and surface-to-volume ratio.

## **Batch Selection**

The protocol covers three primary batches manufactured at commercial scale using the approved process. Batch genealogy, manufacturing dates, and packaging lots are recorded in the batch appendix. All batches were released against the approved specification prior to set-down.

## **Data Evaluation and Reporting**

Stability data are evaluated per ICH Q1E. Regression analysis is applied to quantitative attributes exhibiting change over time; poolability across batches is assessed at the 0.25 significance level. Annual reports summarize all data generated in the reporting interval.

## **Deviations**

Chamber excursion of +3°C for 15 hours affected the 3-month interval for batch 001; mean kinetic temperature remained within limits.

## **Acceptance Criteria**

The stability acceptance criterion applied at each interval is — Any individual unspecified degradant: not more than 0.20%. All other registered specification tests apply unchanged.

## **Data Evaluation and Reporting**

Stability data are evaluated per ICH Q1E. Regression analysis is applied to quantitative attributes exhibiting change over time; poolability across batches is assessed at the 0.25 significance level. Annual reports summarize all data generated in the reporting interval.

## **Photostability**

Confirmatory photostability testing per ICH Q1B was completed during development and is not repeated in this protocol. The product is labeled for storage in the original container to protect from light.

## **Container Closure and Orientation**

Samples are stored in the commercial container closure system. Bottles are stored upright; a subset of inverted samples is included at the initial and final timepoints to assess closure-liner interaction. Induction seal integrity is verified at each pull.