

Gestrinal — Stability Study Design

Protocol STAB-P-2021-18

Photostability

Confirmatory photostability testing per ICH Q1B was completed during development and is not repeated in this protocol. The product is labeled for storage in the original container to protect from light.

Storage Conditions

The primary (long-term ICH condition) storage condition for this protocol is $25^{\circ}\text{C} \pm 2^{\circ}\text{C}$ / $60\% \text{ RH} \pm 5\% \text{ RH}$. Accelerated studies at $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$ / $75\% \text{ RH} \pm 5\% \text{ RH}$ were completed through six months during development and are not part of this protocol.

Pull Schedule

Interval (months)	0	3	6	9	12	18	24
Long-term testing	X	X	X	X	X	X	X

Bracketing Justification

Bracketing is justified under ICH Q1D on the basis that the 10 mg and 100 mg presentations bracket the 40 mg presentation with respect to fill volume and surface-to-volume ratio.

Stability Acceptance Criteria

The stability acceptance criterion applied at each interval is — Assay: 95.0–105.0% of label claim. All other registered specification tests apply unchanged.

Batch Selection

The protocol covers three primary batches manufactured at commercial scale using the approved process. Batch genealogy, manufacturing dates, and packaging lots are recorded in the batch appendix. All batches were released against the approved specification prior to set-down.

Data Evaluation and Reporting

Stability data are evaluated per ICH Q1E. Regression analysis is applied to quantitative attributes exhibiting change over time; poolability across batches is assessed at the 0.25 significance level. Annual reports summarize all data generated in the reporting interval.

Protocol Deviations

Chamber excursion of +3°C for 9 hours affected the 6-month interval for batch 001; mean kinetic temperature remained within limits.

Post-Approval Commitments

As a condition of approval, three consecutive production batches of Gestralin will be entered into the stability program, and the retest period will be re-evaluated after 36-month data are available.

Analytical Procedures

Assay and degradation products are determined by the stability-indicating reversed-phase HPLC procedure described in the analytical methods dossier. Dissolution testing uses USP Apparatus II at 50 rpm. All procedures were validated in accordance with ICH Q2(R1) prior to protocol initiation.

Batch Selection

The protocol covers three primary batches manufactured at commercial scale using the approved process. Batch genealogy, manufacturing dates, and packaging lots are recorded in the batch appendix. All batches were released against the approved specification prior to set-down.

Data Evaluation and Reporting

Stability data are evaluated per ICH Q1E. Regression analysis is applied to quantitative attributes exhibiting change over time; poolability across batches is assessed at the 0.25 significance level. Annual reports summarize all data generated in the reporting interval.