

Fluventra — Stability Protocol

Protocol STAB-P-2019-76

Sample Inventory and Pull Windows

Sample inventory per condition and timepoint is defined in the pull schedule appendix. Scheduled pulls are performed within ± 7 days of the nominal date; pull records are maintained in the site stability management system.

Storage Conditions

The primary (zone IVb condition) storage condition for this protocol is $30^{\circ}\text{C} \pm 2^{\circ}\text{C}$ / $75\% \text{ RH} \pm 5\% \text{ RH}$. Accelerated studies at $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$ / $75\% \text{ RH} \pm 5\% \text{ RH}$ were completed through six months during development and are not part of this protocol.

Testing Schedule

Interval (months)	0	1	3	6	9	12	18
Long-term testing	X	X	X	X	X	X	X

Protocol Deviations

The 9-month samples for batch 003 were held at ambient conditions for 18 hours during chamber maintenance before testing; trending showed no anomalous results.

Acceptance Criteria

The stability acceptance criterion applied at each interval is — Dissolution: $Q = 80\%$ in 30 minutes. All other registered specification tests apply unchanged.

Bracketing and Matrixing

A bracketing design per ICH Q1D tests the extremes of the strength range (10 mg and 100 mg); the 20 mg strength shares identical formulation composition and packaging and is therefore not tested.

Stability Commitments

As a condition of approval, three consecutive production batches of Fluventra will be entered into the stability program, and the retest period will be re-evaluated after 36-month data are available.

Data Evaluation and Reporting

Stability data are evaluated per ICH Q1E. Regression analysis is applied to quantitative attributes exhibiting change over time; poolability across batches is assessed at the 0.25 significance level. Annual reports summarize all data generated in the reporting interval.

Container Closure and Orientation

Samples are stored in the commercial container closure system. Bottles are stored upright; a subset of inverted samples is included at the initial and final timepoints to assess closure-liner interaction. Induction seal integrity is verified at each pull.