

Kolzepine — Post-Approval Stability Protocol

Protocol STAB-P-2020-78

Analytical Procedures

Assay and degradation products are determined by the stability-indicating reversed-phase HPLC procedure described in the analytical methods dossier. Dissolution testing uses USP Apparatus II at 50 rpm. All procedures were validated in accordance with ICH Q2(R1) prior to protocol initiation.

Container Closure and Orientation

Samples are stored in the commercial container closure system. Bottles are stored upright; a subset of inverted samples is included at the initial and final timepoints to assess closure-liner interaction. Induction seal integrity is verified at each pull.

Study Storage Conditions

The primary (refrigerated condition) storage condition for this protocol is $5^{\circ}\text{C} \pm 3^{\circ}\text{C}$. Accelerated studies at $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$ / $75\% \text{ RH} \pm 5\% \text{ RH}$ were completed through six months during development and are not part of this protocol.

Pull Schedule

Long-term samples are tested at 0, 3, 6, 12, 18, 24, 36 and 48 months.

Stability Acceptance Criteria

The stability acceptance criterion applied at each interval is — Assay: 95.0–105.0% of label claim. All other registered specification tests apply unchanged.

Post-Approval Commitments

A post-approval commitment is made to enrol one commercial batch of Kolzepine per year into the ongoing stability program and to report confirmed out-of-specification results within 30 days.

Analytical Procedures

Assay and degradation products are determined by the stability-indicating reversed-phase HPLC procedure described in the analytical methods dossier. Dissolution testing uses USP Apparatus II at 50 rpm. All procedures were validated in accordance with ICH Q2(R1) prior to protocol initiation.

Data Evaluation and Reporting

Stability data are evaluated per ICH Q1E. Regression analysis is applied to quantitative attributes exhibiting change over time; poolability across batches is assessed at the 0.25 significance level. Annual reports summarize all data generated in the reporting interval.

Bracketing and Matrixing

In accordance with ICH Q1D, bracketing is applied: the 20 mg and 40 mg strengths are tested, and the intermediate 25 mg strength is covered because all strengths are direct scale-equivalents of a common blend in the same container closure system.

Deviations

No deviations from this protocol have been recorded to date.

Photostability

Confirmatory photostability testing per ICH Q1B was completed during development and is not repeated in this protocol. The product is labeled for storage in the original container to protect from light.

Batch Selection

The protocol covers three primary batches manufactured at commercial scale using the approved process. Batch genealogy, manufacturing dates, and packaging lots are recorded in the batch appendix. All batches were released against the approved specification prior to set-down.

Container Closure and Orientation

Samples are stored in the commercial container closure system. Bottles are stored upright; a subset of inverted samples is included at the initial and final timepoints to assess closure-liner interaction. Induction seal integrity is verified at each pull.