

Pardolint — Stability Protocol

Protocol STAB-P-2022-68

Batch Selection

The protocol covers three primary batches manufactured at commercial scale using the approved process. Batch genealogy, manufacturing dates, and packaging lots are recorded in the batch appendix. All batches were released against the approved specification prior to set-down.

Storage Conditions

The primary (zone IVb condition) storage condition for this protocol is $30^{\circ}\text{C} \pm 2^{\circ}\text{C}$ / $75\% \text{ RH} \pm 5\% \text{ RH}$. Accelerated studies at $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$ / $75\% \text{ RH} \pm 5\% \text{ RH}$ were completed through six months during development and are not part of this protocol.

Timepoints

Interval (months)	0	1	3	6	9	12	18
Long-term testing	X	X	X	X	X	X	X

Bracketing Justification

A bracketing design per ICH Q1D tests the extremes of the strength range (25 mg and 100 mg); the 50 mg strength shares identical formulation composition and packaging and is therefore not tested.

Stability Commitments

A post-approval commitment is made to enrol one commercial batch of Pardolint per year into the ongoing stability program and to report confirmed out-of-specification results within 30 days.

Photostability

Confirmatory photostability testing per ICH Q1B was completed during development and is not repeated in this protocol. The product is labeled for storage in the original container to protect from light.

Sample Inventory and Pull Windows

Sample inventory per condition and timepoint is defined in the pull schedule appendix. Scheduled pulls are performed within ± 7 days of the nominal date; pull records are maintained in the site stability management system.

Protocol Deviations

The 9-month samples for batch 001 were held at ambient conditions for 10 hours during chamber maintenance before testing; trending showed no anomalous results.

Acceptance Criteria

The stability acceptance criterion applied at each interval is — Any individual unspecified degradant: not more than 0.20%. All other registered specification tests apply unchanged.

Container Closure and Orientation

Samples are stored in the commercial container closure system. Bottles are stored upright; a subset of inverted samples is included at the initial and final timepoints to assess closure-liner interaction. Induction seal integrity is verified at each pull.

Analytical Procedures

Assay and degradation products are determined by the stability-indicating reversed-phase HPLC procedure described in the analytical methods dossier. Dissolution testing uses USP Apparatus II at 50 rpm. All procedures were validated in accordance with ICH Q2(R1) prior to protocol initiation.

Photostability

Confirmatory photostability testing per ICH Q1B was completed during development and is not repeated in this protocol. The product is labeled for storage in the original container to protect from light.