

Velmorapt — Commercial Stability Protocol

Protocol STAB-P-2024-70

Sample Inventory and Pull Windows

Sample inventory per condition and timepoint is defined in the pull schedule appendix. Scheduled pulls are performed within ± 7 days of the nominal date; pull records are maintained in the site stability management system.

Study Storage Conditions

The primary (long-term ICH condition) storage condition for this protocol is $25^{\circ}\text{C} \pm 2^{\circ}\text{C}$ / $60\% \text{ RH} \pm 5\% \text{ RH}$. Accelerated studies at $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$ / $75\% \text{ RH} \pm 5\% \text{ RH}$ were completed through six months during development and are not part of this protocol.

Pull Schedule

Interval (months)	0	3	6	9	12	18	24
Long-term testing	X	X	X	X	X	X	X

Bracketing Justification

Bracketing is justified under ICH Q1D on the basis that the 5 mg and 40 mg presentations bracket the 25 mg presentation with respect to fill volume and surface-to-volume ratio.

Post-Approval Commitments

The first three commercial-scale batches of Velmorapt will be placed on the long-term stability program, with at least one production batch added annually thereafter.

Deviations

No deviations from this protocol have been recorded to date.

Stability Acceptance Criteria

The stability acceptance criterion applied at each interval is — Assay: 95.0–105.0% of label claim. All other registered specification tests apply unchanged.

Batch Selection

The protocol covers three primary batches manufactured at commercial scale using the approved process. Batch genealogy, manufacturing dates, and packaging lots are recorded in the batch appendix. All batches were released against the approved specification prior to set-down.

Analytical Procedures

Assay and degradation products are determined by the stability-indicating reversed-phase HPLC procedure described in the analytical methods dossier. Dissolution testing uses USP Apparatus II at 50 rpm. All procedures were validated in accordance with ICH Q2(R1) prior to protocol initiation.